


# MAX-606



## PATIENT POSITIONING LASER SYSTEM



CEMARLINE



# MAXX-606

## MAXX-606 Laser Specifications

- Four Manual Axis Movements:
  - X Axis (Horizontal Movement): 20mm (0.8")  
+10mm / -10mm (from Isocenter)
  - Y Axis (Vertical Movement): 20mm (0.8")  
+10mm / -10mm (from Isocenter)
  - Roll Axis (Horizontal Rotation): 10° +5° / -5°
  - Pitch Axis (Vertical Rotation): 6° +3° / -3°
- Crosshair Beam Rotation: 15° +7.5° / -7.5°
- Focus Adjustment: 0.5m - 3.0m (10ft)
- Crosshair Beam Fan Out: 60° Both Axis
- Available Colors: **Green** or **Red**
- Laser Power Adjustment: < 1mW
- Maximum Line Width @ 3m (10'): < 1.2mm (0.047")
- Operating Temperature: 15°C - 35°C (10°F - 95°F)
- Power Requirements: 110V AC – 240V AC, 0.8A, 50Hz - 60Hz
- Weight: 1.8kg (3.9 lbs.)
- Dimensions:
  - Length: 210mm (8.2")
  - Width: 118mm (4.6")
  - Height: 100mm (3.9")



### Cemar Electro Inc. Headquarters and International Sales

Location: 1370 55e Avenue, Lachine, Quebec, H8T 3J8, Canada  
Email Address: [cemar@cemarelectro.com](mailto:cemar@cemarelectro.com)  
Telephone Number: 1-514-831-5807

### Cemar Medical Laser SAS EU Authorized Representative

Location: 1 Avenue, Christian Doppler, 77700, France  
Email Address: [philippe.zimmer@cemar-medical-laser.com](mailto:philippe.zimmer@cemar-medical-laser.com)  
Telephone Number: 33-1-8110-5172

### MIS Healthcare UK Authorized Representative

Location: 26-27 Capitol Park, London, NW9 0EQ, United Kingdom  
Email Address: [abritton@mishealthcare.co.uk](mailto:abritton@mishealthcare.co.uk)  
Telephone Number: 44-208-205-9500

