


MAXX-606



PATIENT POSITIONING LASER SYSTEM



—CEMARLINE— 

MAXX-606

MAXX 606 Laser Specifications

- **Four Manual Axis Movements:**
 - X Axis (Horizontal Movement): 20mm (0.8") $\pm$ 10mm (from Isocenter)
 - Y Axis (Vertical Movement): 20mm (0.8") $\pm$ 10mm (from Isocenter)
 - Roll Axis (Horizontal Rotation): $10^{\circ} \pm 5^{\circ}$
 - Pitch Axis (Vertical Rotation): $6^{\circ} \pm 3^{\circ}$
- Crosshair Beam Rotation: $15^{\circ} \pm 7.5^{\circ}$
- Focus Adjustment: 0.5m - 3.0m (1.6ft - 10ft)
- Crosshair Beam Fan Out: 60° Both Axis
- Available Colors: **Green** or **Red**
- Laser Power Adjustment: <1 mW
- Maximum Line Width at 3m (10'): <1.2 mm (0.047")
- Operating Temperature: $5^{\circ}\text{C} - 35^{\circ}\text{C}$ ($41^{\circ}\text{F} - 95^{\circ}\text{F}$)
- Power Requirements: 110V AC – 240V AC, 0.8A, 50Hz – 60Hz
- **Dimensions:**
 - Length: 210mm (8.2")
 - Width: 118mm (4.6")
 - Height: 100mm (3.9")
 - Weight: 1.8 kg (3.9 lbs)



Cemar Electro Inc.
Head Office and International Sales
1370 55^e Avenue, Lachine, Quebec,
H8T 3J8, Canada
Email Address: cemar@cemarelectro.com
Telephone Number: +1-514-631-5807

Cemar Medical Laser SAS
European Representative - CE
14 Avenue de L'Europe,
77144 Montevrain, France
Email Address: philippe.zimmer@cemar-medical-laser.com
Telephone Number: +336 0778 5626

WWW.CEMARELECTRO.COM

Doc: 606-08-26

CEMARLINE

