

MAXX-700

PATIENT POSITIONING LASER SYSTEM



CEMARLINE





MAXX-700



MAXX-700 Laser System Features

- Two Auto Alignment Verification Sensors
- Auto Ambient Light Sensor Adjustment
- Four Automatic Axis Movements
- Macro and Micro Adjustments
- Motor Drives with Precision
- Available in **Green** or **Red**

Optional Mounting Brackets



CEMARLINE




MAXX-700

MAXX Tablet Operations




- Communicates with MAXX-700, MAXX-1100, and MAXX Combo Set Configurations
- LCD Touch Screen Tablet with Protective Case
- Durable Wireless Communication
- Updated User-Friendly Interface
- Extended Battery Capacity

MAXX-700 Application Start Page



MAXX-700 Laser Position Page



MAXX-700 Laser Settings Page

CEMARLINE

MAXX-700

MAXX-700 Laser Specifications

- **Four Automatic Axis Movements:**
 - o X Axis (Horizontal Movement): 20 mm (0.8") $\pm$ 10 mm (from Isocenter)
 - o Y Axis (Vertical Movement): 20 mm (0.8") $\pm$ 10 mm (from Isocenter)
 - o Roll Axis (Horizontal Rotation): 10°, $\pm$ 5°
 - o Pitch Axis (Vertical Rotation): 6° $\pm$ 3°
- Movement and Position Detection
- Crosshair Beam Rotation: 15° $\pm$ 7.5°
- Focus Adjustment: 0.5 m - 3.0 m (1.6t - 10ft)
- Crosshair Beam Fan-Out: 60° on Both Axes
- Precision Motor-Driven Displacements
- Available Colors: **Green** or **Red**
- Laser Power Adjustment: < 1 mW
- Maximum Line Width at 3 m (10'): < 1.0 mm (0.04")
- The Tablet Can Communicate with up to Five Lasers per Room
- Operating Temperature: 5°C - 35°C (41°F - 95°F)
- Power Requirements: 100 V AC – 240 V AC, 0.8 A, 50 Hz – 60 Hz
- **Dimensions:**
 - o Length: 260 mm (10.24")
 - o Width: 116 mm (4.5")
 - o Height: 94 mm (3.7")
 - o Weight: 2.9 kg (6.4 lbs.)



Cemar Electro Inc.

Headquarters and International Sales

Location: 1370 55^e Avenue, Lachine, Quebec,
H8T 3J8, Canada

Email Address: cemar@cemarelectro.com

Telephone Number: 1-514-631-5807

WWW.CEMARELECTRO.COM

Cemar Medical Laser SAS

European representative

14 Avenue de L'Europe,
77144 Montevrain, France

Email Address: philippe.zimmer@cemar-medical-laser.com

Telephone Number: +336 0778 5626

Doc: 700-08-26

CEMARLINE

